



www.technosource.in

INSTRUCTION MANUAL OF DNA/RNA Isolation KIT

(KT-NA001/ KT-NA005/ KT-NA010)

Manufactured and Marketed by:

TechResource

188, Nagdevi Street, 2nd Floor, Mumbai 400003.

Tel: +91 9819095429 Email: sales@techresource.in

Kit components

Sr.no.	Components	Store at
1.	Buffered Agarose (1%)	4°C
2.	5X Tank buffer	4°C
3.	Loading cum dissolution dye with EtBr	4°C
4.	EtBr (10 mg/mL)	4°C
5.	KMnO ₄ (Powder)	RT
6.	D/W	RT
7.	Homogenization buffer	RT
8.	Sodium acetate (Precipitation buffer)	RT
9.	TE buffer	RT
10.	Chaotropic buffer	RT
11.	Chloroform	4°C
12.	Ethanol	RT
13.	70% Ethanol	RT

Equipments and consumables required

but not provided in the kit:

- Fresh Chicken brain (1g for 1 Expt.)
- Mortar & pestle
- Nylon filter pouch
- Ice bath
- Hot water bath
- Refrigerator
- Beakers (5, 15 mL): 1 no. each
- Centrifuge (for microfuge tube)
- Centrifuge tubes
- Horizontal gel electrophoresis apparatus [maximum gel size 7 (w) X 10 (l) cm]
- One 100 mL measuring cylinder
- Two 250 mL beakers
- Power supply
- UV Transilluminator
- Micropipettes-
 - 0.5-10 µL
 - 5-50 µL
 - 10-100 µL
 - 100-1000 µL
- Micropipette Tips

Horizontal electrophoresis apparatus

Electrophoresis apparatus manufactured those marked as * are suitable for this kit.

- 1) ***Sleeksub** (TechnoSource) Gel size - 7 (w) X 10 (l) Cat. no. 08-00102/105
- 2) **Gelmariner** (TechnoSource) Gel size - 10 (w) X 15 (l) Cat. No. 10-00152/155/156
- 3) **Greyhound** (TechnoSource) Gel size - 15 (w) X 20 (l) Cat. no. 09-00202/5/6/3
- 4) **Submariner** (TechnoSource) Gel size - 23 (w) X 20 (l) Cat. no. 17-00232/5/6/3
- 5) **Greywhale** (TechnoSource) Gel size - 30 (w) X 30 (l) Cat. no. 12-00302/5/6/3
- 6) ***Slimsub** (EcoSep) Gel size - 3 (w) X 10 (l) Cat. no. E2100031/32
- 7) ***Ecosub** (EcoSep) Gel size - 7 (w) X 10 (l) Cat. no. E2200071/72
- 8) **Midisub** (EcoSep) Gel size - 10 (w) X 10 (l) Cat. no. E2400101/102
- 9) **Tubbysub** (EcoSep) Gel size - 15 (w) X 10 (l) Cat. no. E2300151/152

Preparation instructions:

0.5X tank buffer:

0.5X tank buffer is to be prepared from 5X tank buffer by diluting with D/W. Refer to the table below

Table 1:

	1 Expt.	5 Expts.	10 Expts.
D/W	270 mL	1350 mL	2700 mL
5X Tank buffer	30 mL	150 mL	300 mL
Total	300 mL	1500 mL	3000 mL

Store 0.5X tank buffer at 4°C.

Choice between TAE and TBE

The choice between TAE and TBE can be made based on the following properties

- a) TAE has low buffering capacity and will get exhausted if electrophoresis is carried out for prolonged times.
- b) TBE is more expensive than TAE.
- c) Double stranded DNA fragments migrate faster through TAE than through TBE.
- d) The resolving power of TAE is better for higher molecular weight and for supercoiled DNA. TBE is better for lower-molecular weight DNA.

Choice of Agarose concentrations

Agarose conc in gel (% w/v)	Range of linear DNA	V grad-ient (V/cm)	Notes
0.3	5-60	5	Very fragile. Opt for PAGE if possible
0.6	1-20	5	
0.7	0.8-10	5	BPB runs at approx. 1 kb, XC at approx. 4 kb
0.8	0.5-7	10	
0.9	0.5-7	10	BPB runs at approx. 0.5 kb
1.0	0.5-6	10	
1.2	0.4-6	10	Pour agarose at 55°C in pre warmed tray
1.5	0.2-3	10	
2.0	0.1-2	20	
3.0	0.1-1	20	

Troubleshooting Guide 1

Indication	Cause	Remedy
No bands/very faint bands of DNA	Less tissue taken	Use more tissue for extraction
	Incomplete homogenization	Homogenize tissue thoroughly in homogenization buffer
	Heating not enough or overheated the homogenized tissue	Heat the homogenized tissue at proper temperature as mentioned in the protocol
	Poor EtBr staining	Use proper concentration of EtBr (i.e. 0.5 µg/mL)
Bands are fuzzy	Too much tissue was taken	Weigh 1 g of tissue as mentioned in the protocol
	DNA/RNA contaminated with proteins	Add appropriate amount of chloroform and spun vigorously
No bands/very faint bands of RNA	Nuclease contamination	Make sure all the glassware, plastic-ware and solutions are autoclave before use
	Incomplete homogenization	Homogenize tissue thoroughly in homogenization buffer
	Poor EtBr staining	Use proper concentration of EtBr (i.e. 0.5 µg/mL)
Incomplete ethanol evaporation	Allow the ethanol to evaporate completely from RNA pellet before dissolving in TE buffer	

Troubleshooting Guide 2

Indication	Cause	Remedy
1. Problems with power supply		
Power supply will not switch 'ON', i.e. light on mains switch/LED indicator does not come 'ON'.	Check 230 V main supply	Restore the mains supply
	Power supply fuse blown	Try replacing with correct fuse only once
	Fuse blows again	Contact the manufacturer of the power supply
2. Problems with horizontal electrophoresis unit		
Power supply 'ON' but no movement of samples/ bubbles in the buffer	Cathode and anode are not connected to the right output terminals.	Connect the cathode and anode to the right output terminals i.e. red to red and black to black
	Electrical connecting cords are not properly inserted.	Insert the electrical connecting cords properly.
	The connecting cords are damaged.	Use new connecting cords
	Low level of buffer in the apparatus	Raise the level of buffer in the apparatus such that the circuit is formed.
No bands seen	Presence of cello tape in tray	Remove cello tape and check
	Oxidation/deposit formation on the contact terminals/ connectors	Disconnect the instrument. Clean the terminals. Reconnect the instrument and try again.
DNA was electrophoresed off the gel	DNA was electrophoresed off the gel	Electrophorese the gel for less time, at lower voltage
	To low concentration of sample	Concentrate the sample before loading
	sample was degraded	Avoid nuclease contamination
	For ethidium bromide-stained DNA, improper UV light source was used.	Use short-wavelength (254 nm) UV light for greater sensitivity.

DNA/RNA isolation and Separation of DNA/RNA by Agarose Gel Electrophoresis

Introduction:

The isolation of DNA from any source usually begins with lysis, or breakdown of tissue or cells. Lysis is carried out in the salt solution containing detergent such as SDS to denature proteins. Lysis of the tissue is a crucial step which allows nucleic acids to be released from the nucleus by disintegrating it. Heating (at 50-60°C) during lysis step enhances the process and results in breakdown of cells. Once the lysis is complete, the lysate is treated with chloroform to remove all the proteins and lipids from the lysate. DNA is recovered from the aqueous solution by treating with ethanol. The isolation of RNA involves four basic steps. First, lysis of the cells by addition of lysis solution containing detergent. After cell has lysed, treat the lysate with chloroform to separate

RNA from the solution. RNA is recovered from the aqueous solution by addition of cold ethanol. Precipitated RNA is then washed with 70% ethanol to remove any salt or impurities from the solution.

Theory:

Electrophoresis is a technique in which charged molecules migrate from cathode to anode under influence of electric current. Agarose is a polysaccharide obtained from red algae. Agarose forms inert matrix through which charged molecules migrate towards anode. Smaller molecules move faster compared to larger ones, in a process called sieving. Higher the percentage of agarose lower is the pore size and vice versa. Most agarose gels are made between 0.7% and 2% of agarose depending on size of DNA/RNA to be separated. The DNA/RNA can be visualized by adding ethidium bromide (EtBr) in the gel or loading buffer.

DNA isolation protocol

- 1) Sterilize all the glassware by autoclaving.
- 2) Weigh 1 g of fresh chicken brain tissue.
- 3) Add 1.25 mL homogenization buffer and 190 μ L Chaotropic buffer. Homogenize the chicken using mortar and pestle.
- 4) Transfer the paste in a small (5 mL) beaker. Heat at 50-60°C in water bath for 10 mins.
- 5) Add equal volume of chloroform (i.e. 1.5 mL) and mix well.
- 6) Transfer the content in a centrifuge tube and spin hard for 10-15 mins.
- 7) Collect supernatant in a chilled beaker (15 mL) or tube and cool in freezer/keep on ice bath for 10 mins.
- 8) Add 2 ½ volume of ethanol (i.e if volume is 3 mL then add 7.5 mL of ethanol) to the supernatant from sides of the beaker/tube very slowly.
- 9) Swirl gently and cool in freezer/keep on ice bath for 10 mins.
- 10) Collect DNA fibers in a microfuge tube with a glass rod. Again swirl the content gently and keep in freezer/ ice bath. Collect the DNA fibers in a same microfuge tube.
- 11) Add 0.5 mL TE buffer and store at 4°C.

RNA isolation protocol

- 1) Sterilize all the glassware by autoclaving.
- 2) Weigh 1 g of fresh chicken brain tissue and place tissue in nylon filter pouch.
- 3) Immerse the bag of tissue in a mortar containing 2 mL homogenization buffer.
- 4) Mash the tissue by crushing it within bag. Empty bag contents in mortar.
- 5) Add 300 μ L precipitation buffer, mix and add equal volume of chloroform. Mix vigorously.

6) Transfer the content in sterile microfuge tube and spin hard for 15 mins.

7) Remove the supernatant in sterile microfuge tube. Again spin to remove debris and fibres.

8) Collect the supernatant in tube and add 2.5 mL cold ethanol. Mix thoroughly and chill for 10 mins. on ice bath.

9) Spin to collect the precipitate.

10) Wash the precipitate with 70% ethanol. Spin for 15 mins.

11) Discard the supernatant and keep the tube open at RT until the ethanol has evaporated.

12) Re-suspend the precipitate in 0.5 mL TE. Store at 4°C.

Electrophoresis protocol

1) Rinse all items, beakers, measuring cylinders, pipettes and electrophoretic apparatus with distilled water, dry and use.

2) Prepare 0.5 X tank buffer as mentioned in table 1. Cover with a plastic sheet and chill in the refrigerator.

Sealing the tray:

3) Tape both the ends of the tray making sure that the tape is high enough and place it on a levelled surface.

Casting the gel:

4) If using the comb stand: place the comb stand on the tray/ If using the fixed height comb: insert the comb in the slots provided.

5) Melt the agarose in a hot water bath.

6) Cool the agarose to approximately 50 °C and add EtBr (final concentration 0.5 μ g/ml) if desired. Pour the warm solution into the tray to a gel thickness of about 4-5 mm. Allow the gel to set for minimum 15 mins.

7) Carefully remove the tape and place the gel tray in the apparatus such that the wells are closer to the cathode (black connector).

8) Pour 0.5X tank buffer into the apparatus such that buffer should be 2-4 mm above the gel.

9) Remove the comb gently.

Sample preparation:

10) Mix sample and loading dye in a proportion of 1:5 (i.e. for 5 μ L of sample add 1 μ L of dye)

Running the gel:

11) Load 20 μ L of the sample per well.

12) Connect the apparatus to power supply. Start the power supply. Adjust the voltage and current as per your requirement. (For a small apparatus of 3 x 10 cm adjust to a constant voltage of 50-70) *Check if the samples are migrating in the correct direction.*

13) Switch the power supply off when the blue tracking dye reaches towards the end of the gel.

14) Remove the gel tray from the apparatus, drain out all the buffer carefully and dry the bottom surface of the gel tray with blotting paper. If the tray is UV transparent, place the gel tray directly over the window of trans-illuminator. If the tray is UV opaque, then remove the gel from the tray and place over the window of trans-illuminator/ gel documentation system.

15) Switch on the trans-illuminator and visualize the bands.

16) Discard all buffers containing the nucleic acid intercalating dye into a beaker containing purple crystals of dye oxidiser.